

☎ (0341) 3021374

🌐 info@gbi-medika.com

📍 Jl. Band. Abd. Rahman Saleh No.GG
No. 3, Cemorokandang, Kec.
Kedungkandang, Kota Malang, Jawa
Timur 65138

📷 @gbimedikbio



MEDIK BIO® From Research To Innovation

www.medikbioindonesia.com

MEDIK BIO®
From Research To Innovation

Distributed by:
PT GAMA BIOTEK INDONESIA





INNOVATION OF MEDIKBIO

Design Innovation



MEDIKBIO medical designer have always maintain close communication with medical service providers, to help them providing better solutions to their patents. With long-standing experience in product structure and functional experience, MEDIKBIO continue to forget ahead with innovation and play a vital role in helping medical service providers to improve the quality of care

Research Innovation



MEDIKBIO medical owns independent research an development center. since its inception every year MEDIKBIO invest 10% of its annual sales in R&D, launched 10 new product, and applied at least 10 patented technologies in an average basis, which lead MEDIKBIO to continuously define the quality if medical product and new standards of efficiency

Manufacturing Innovation



MEDIKBIO medical, as a pioneer in the field of medicaal furniture, always choose branded suppliers strctly to ensure the the safety and reliability of raw material. MEDIKBIO adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operaator, rigorous quality control, only to guarantee MEDIKBIO's commitment on quality

Maintenance Innovation



MEDIKBIO medical have more than 50 engineers, and qualified service term that composed of professional subcontractors, trained technicians, and after sale staff. They are experienced and have excellent skills, not only to respond each service request on time, but also precaution every detail in equipment operation, always focus on product safely

Service Innovation



MEDIKBIO medical have thorough service architecture: 7 days x 24hours online customer call canter, complate spare parts warehouse, advance CRM system, precautionary service, etc. No matter where customer is, MEDIKBIO service will be always close at hand



Doctor-patient communication system

Rehabilitation guidance, follow-up, tracking, reminders, allowing patients to communicate with doctors through text, voice, and video anytime, anywhere, creating a mobile medical platform for in-depth doctor-patient communication.



Real-time vital sign monitoring

Monitor heart rate, blood pressure, blood oxygen, respiration, temperature, pulse, etc.



Remote Monitoring System

Real-time video monitoring and audio-video conversation.




**Vital Sign Monitoring
Patient 360 Monitoring
Infusion Monitoring
Video Monitoring**




Equipped with Smart Hospital Bed System

Clinical data collection, screening, integration. Generate medical big data to assist decision-making.



Control system





75°±10°
Back-rest

35°±10°
Knee-rest

400mm
Height
Adjustment
530-930mm

15°±2°
Reverse
Trendelenburg

15°±2°
Trendelenburg

17°±2°
Lateral Tilt
Left

17°±2°
Lateral Tilt
Right

weight
Weighing System

CPR
Instant Release

Cardiac
Chair Position

0-90°
Angle Indicator

75+15mm
Auto-regression

30°
Position Line

LINAK
Motor

Brake Indicator

Rail Lock
Indicator

Z7z
ICU Smart Bed
SKZH-ABCZ

Technical Parameters

External Length	2127mm
External Width	1070mm
Bed Height Min.	530mm
Bed Height Max.	930mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

✓ LINAK Motor	2pcs
✓ Telescopic Column Motor	3pcs
✓ 6" Double Side Secure® Castors, Central Locking	4pcs
✓ PP Side Rail with Integrated LCD Control Panel	1set
✓ PP Bed Ends with Safe Lock	1set
✓ X-ray Transparent Phenolic Resin Bed Platform	1set
✓ Patent Angle Indicator	4pcs
✓ Φ85mm Wall Bumper Wheel	4pcs
✓ Trapeze Socket	2pcs
✓ IV Pole	1set
✓ IV Pole Prevision	8pcs
✓ Drainage Hook	4pcs
✓ Central Locking Pedal	2pcs
✓ PP Mattress Retainer	6pcs
✓ Nurse Cotroller	1pc
✓ Emergency Stop Button	1pc
✓ Night Light	1set
✓ Manual CPR	
✓ Electric CPR	
✓ Weighing System	



70°±10°
Back-rest

28°±10°
Knee-rest

450mm
Height
Adjustment
400-850mm

17°±2°
Reverse
Trendelenburg

17°±2°
Trendelenburg

Cardiac
Chair Position

200mm
Bed Extension

weight
Weighing System

CPR
Instant Release

0-90°
Angle Indicator

67±57mm
Auto-regression

30°
Position Line

Built-in LED
Torch

LINAK
Motor

Nurse Call

Siderail
Unlocked

Brake
Unlocked

Lowest
Position

Out-of-bed
Position Warnings

C8c
Electric ICU Bed
SKD-CV

Technical Parameters

External Length	2180mm
External Width	995mm
Bed Height Min.	400mm
Bed Height Max.	850mm
Mattress Platform	1995×860mm
Safe Working Load	250kg

Technical Configuration

✓ LINAK Motor	4pcs	✓ Nurse Cotroller	1pc
✓ 5" Double Side Castors with directional wheel function	4pcs	✓ Cardiac Chair Position	
✓ PP Removable C-series Headboard and Footboard Φ113 Bumper Wheels	1set	✓ Electric CPR	
✓ PP 4-pcs C-Series Siderails	1set	✓ Manual CPR	
✓ Head-end Siderails with Integrated Controller and Warning Lights	1set	✓ Auto Break at 30°	
✓ Foot-end Siderails with Glasses Holder, LCD Screen and 15W Wireless Charger	1set	✓ Emergency Stop Button	
✓ PP Detachable Bed Platform With Integrated Mattress Retainer	1set	✓ Night Light	
✓ Mechanical Lever on Leg Board for Calf Raise	1set	✓ Bed Extension	
✓ Backrest with X-ray Cassette Holder	1pc	✓ Weighing System	
✓ Patent Angle Indicator for Back-rest and Trendelenburg Position	4pcs	✓ Bed End with Projection Lights	
✓ IV Pole 4 Hooks	1set	✓ Nurse Call	
✓ IV Pole Prevision	4pcs	✓ Siderail Unlocked	
✓ Drainage Hook	2pcs	✓ Brake Unlocked	
✓ Φ32mm Monkey Pole Holes	2pcs	✓ Lowest Position	
✓ Hand Cotroller	1pc	✓ Out-of-bed Position Warnings	



75°±10°
Back-rest

20°±10°
Knee-rest

400mm
Height
Adjustment
550-950mm

18°±2°
Reverse
Trendelenburg

18°±2°
Trendelenburg

weight
Weighing System

CPR
Instant Release

0-90°
Angle Indicator

75mm
Auto-regression

200mm
Bed Extension

Cardiac
Chair Position

LINAK
Motor

X9x
Electric ICU Bed
SKD-CV

Technical Parameters

External Length	2120-2320mm
External Width	1010mm
Bed Height Min.	550mm
Bed Height Max.	950mm
Mattress Platform	1890-2090×820mm
Safe Working Load	250kg

Optional

☐ Monkey Pole

Technical Configuration

✓ LINAK Motor	6pcs
✓ 6" Double Side Castors, Central Locking	4pcs
✓ PP Side Rail with Integrated Control Panel	1set
✓ PP Bed Ends One Step Removable	1set
✓ Ø85mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	4pcs
✓ 5th Directional Wheel Pedal	2pcs
✓ Central Locking Pedal	2pcs
✓ X-ray Translucent Mattress Platform	1set
✓ Back-up Battery	1set
✓ Night Light	1set
✓ Bed Push Handle for Emergency Transfer	1set
✓ Manual CPR	
✓ Electric CPR	
✓ Weighing System With LCD Screen	



Oasis Series

Y8y
Electric ICU Bed
SKD-CV

Technical Parameters

External Length	2150mm
External Width	1060mm
Bed Height Min.	390mm
Bed Height Max.	840mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Optional


- ☐ Monkey Pole
- ☐ Oxygen Cylinder Holder
- ☐ Basket

Advantages

- ★ Double Auto-regression
- ★ Super-low Bed Position 390mm
- ★ PP Bed Platform, Removable & Washable
- ★ Inserted-type PP Head Foot Board, One Step Release With Safe Lock, Patent Design
- ★ PP Safety Guardrail With Two Steps Unlocking System, Patent Design
- ★ Back Rest 30degree Position, Reduce The Chance Of Lung Infection
- ★ Built-in LED Torch

Technical Configuration

✓ LINAK Motor	4pcs
✓ 5" Double Side Castors, Central Locking System	4pcs
✓ PP Oasis Side Rail Extra Height	1set
✓ PP Oasis Bed Ends with Safe Lock	1set
✓ PP Bed Platform, Removable & Washable	5pcs
✓ X-ray Transparent Phenolic Resin Back Platform	1set
✓ Plastic Socket Inside PP Guardrail for Phone or Glasses	2pcs
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
✓ Ø113mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	6pcs
✓ Three-grade Central Locking Pedal	2pcs
✓ ABS Mattress Retainer	2pcs
✓ Manual CPR	
✓ Electric CPR	
✓ Nurse Controller	
✓ Hand Controller	
✓ Weighing System	
✓ Night Light	
✓ Bed Extension 200mm	
✓ Trapeze Socket	
✓ Back Up Battery	
✓ Fif th and Six th Siderails	



Oasis Series

Y8t
Electric ICU Bed
SKD-C

Technical Parameters

External Length	2140mm
External Width	1060mm
Bed Height Min.	480mm
Bed Height Max.	770mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration



✓ LINAK Motor	4pcs
✓ 5" Double Side Castors, Central Locking System	4pcs
✓ PP Oasis Side Rail Extra Height	1set
✓ PP Oasis Bed Ends with Safe Lock	1set
✓ Stamping Forming Bed Platform	1pc
✓ Plastic Socket Inside PP Guardrail for Phone or Glasses	2pcs
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
✓ Ø113mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs
✓ Central Locking Pedal	2pcs
✓ Emergency Stop Button	1pc
✓ Manual CPR	
✓ Electric CPR	
✓ Lowest position indicator	

Optional

☐ Hand Controller
☐ Weighing System
☐ Nurse Control Panel
☐ Cardiac Chair Position
☐ Bed Extension 20cm
☐ Back-up Battery
☐ Directional Wheels Function
☐ ABS Mattress Retainer
☐ Monkey Pole
☐ Balkan Frame
☐ Linen Holder



360° View



Grace Care
D8d-CPR
Electric Bed
SKD-C

Technical Parameters


External Length	2280mm
External Width	1090mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Optional

- ☐ Monkey Pole
- ☐ Oxygen Cylinder Holder
- ☐ Basket

Technical Configuration

☒ LINAK Motor	4pcs
☒ 5" Double Side Castors, Central Locking	4pcs
☒ PP Grace Side Rail with Integrated Control Panel	1 set
☒ PP Grace Bed Ends with Safe Lock	1 set
☒ Intergrated Control Panel (3.5 inch LCD)	1 pc
☒ Phenolic Resin Backrest platform	1 pc
☒ Trapeze Socket	2pcs
☒ IV Pole 4 Hooks	1 set
☒ IV Pole Prevision	4pcs
☒ Wall Bumper Wheel Ø85mm	4pcs
☒ Three-grade Central Locking Pedal	2pcs
☒ Drainage Hook	2pcs
☒ Bed extension 200mm	
☒ Manual CPR	
☒ Electric CPR	
☒ Emergency Stop Button	
☒ Cassette Holder 448×400mm	
☒ Night Light	
☒ Back-up Battery	



V8v
Electric Bed
SKD-CV

Technical Parameters

External Length	2085mm
External Width	1025mm
Bed Height Min.	420mm
Bed Height Max.	870mm
Mattress Platform	1925×900mm
Safe Working Load	250kg


Optional

- ☐ Monkey Pole
- ☐ Oxygen Cylinder Holder
- ☐ Basket
- ☐ Mechanical Lever

Technical Configuration

✓ LINAK Motor	4pcs	✓ Nurse Cotroller	1pc
✓ 6" Double Side Castors, Central Locking	4pcs	✓ Trapeze Socket	2pcs
✓ PP Side Rail Extra Height	1 set	✓ Cardiac Chair Position	
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1 set	✓ Electric CPR	
✓ PP Bed Platform with Integrated Mattress Retainer	1 set	✓ Manual CPR	
✓ X-ray Translucent Mattress Platform (Back Section)	1 pc	✓ Auto Break at 30°	
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs	✓ Emergency Stop Button	
✓ IV Pole 4 Hooks	1 set	✓ Night Light	
✓ IV Pole Prevision	4pcs	✓ Cassette Holder (Backrest)	
✓ Drainage Hook	2pcs	✓ Bed Extension	
✓ Central Locking Pedal	2pcs	✓ Weighing system	
✓ Hand Cotroller	1 pc		

V6v
Electric Bed
SKD-D



Technical Parameters

External Length	2085mm
External Width	1025mm
Bed Height Min.	430mm
Bed Height Max.	740mm
Mattress Platform	1925×900mm
Safe Working Load	250kg


Technical Configuration

✓ Linear Motor	3pcs
✓ Hand Controller	1pc
✓ 5" Hardshell Castors	4pcs
✓ PP Side Rail Extra Height	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ Patent Angle Indicator For Back-rest Position	2pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

☐ Monkey Pole
☐ Trapeze Socket
☐ Oxygen Cylinder Holder
☐ Basket

V8w
Electric Bed
SKD-C



Technical Parameters

External Length	2085mm
External Width	1020mm
Bed Height Min.	455mm
Bed Height Max.	745mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

✓ LINAK Motor	4pcs
✓ Hand Controller	1pc
✓ 5" Hardshell Castors	4pcs
✓ Foldable Side Rail with ProErgo-Handle™	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

☐ Monkey Pole
☐ Trapeze Socket
☐ Basket

V4w
Manual Bed
SK-CF1



Technical Parameters

External Length	2130mm
External Width	1020mm
Bed Height Min.	475mm
Bed Height Max.	725mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

✓ Manual Cranks	4pcs
✓ 5" Hardshell Castors	4pcs
✓ Foldable Side Rail with ProErgo-Handle™	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

☐ Monkey Pole
☐ Trapeze Socket
☐ Basket

V3v
Manual Bed
SK-CD1



Technical Parameters

External Length	2130mm
External Width	1025mm
Bed Height Min.	430mm
Bed Height Max.	740mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

✓ Manual Cranks	3pcs
✓ 5" Hardshell Castors	4pcs
✓ PP Side Rail Extra Height	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ Patent Angle Indicator For Back-rest Position	2pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

☐ Utility Shelf
☐ Trapeze Bracket/Pole

V2k
Manual Bed
SK-C1




Technical Parameters

External Length	2130mm
External Width	995mm
Bed Height	500mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

☒ Manual Crank	2pcs	☒ IV Pole Prevision	4pcs
☒ 5" Covered Castors	4pcs	☒ Drainage Hook	2pcs
☒ Foldable Side Rail	1 set	☐ Monkey Pole	
☒ PP Bed Ends with Safe Lock with Φ 65 Bumper Wheels	1 set	☐ Trapeze Socket	
☒ PP Bed Platform with Integrated Mattress Retainer	1 set	☐ Basket	
☒ IV Pole 4 Hooks	1 set		

V6k
Electric Bed
SKD-D




Technical Parameters

External Length	2085mm
External Width	995mm
Bed Height Min.	420mm
Bed Height Max.	730mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

☒ Linear Motor	3pcs	☒ IV Pole 4 Hooks	1 set
☒ Hand Controller	1 pc	☒ IV Pole Prevision	4pcs
☒ 5" Covered Castors	4pcs	☒ Drainage Hook	2pcs
☒ Foldable Side Rail	1 set	☐ Monkey Pole	
☒ PP Bed Ends with Safe Lock with Φ 65 Bumper Wheels	1 set	☐ Trapeze Socket	
☒ PP Bed Platform with Integrated Mattress Retainer	1 set	☐ Basket	

V3k
Manual Bed
SK-CD1



Technical Parameters

External Length	2130mm
External Width	995mm
Bed Height Min.	420mm
Bed Height Max.	730mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

☒ Manual Crank	3pcs	☒ IV Pole Prevision	4pcs
☒ 5" Covered Castors	4pcs	☒ Drainage Hook	2pcs
☒ Foldable Side Rail	1 set	☐ Monkey Pole	
☒ PP Bed Ends with Safe Lock with Φ 65 Bumper Wheels	1 set	☐ Trapeze Socket	
☒ PP Bed Platform with Integrated Mattress Retainer	1 set	☐ Basket	
☒ IV Pole 4 Hooks	1 set		

K6k

Electric Bed SKD-D



Technical Parameters


External Length	2105mm
External Width	980mm
Bed Height Min.	400mm
Bed Height Max.	710mm
Mattress Platform	1925×820mm
Safe Working Load	200kg

Technical Configuration

✓ Linear Motor	3pcs	✓ IV Pole Prevision	4pcs
✓ Hand Controller	1pc	✓ Drainage Hook	4pcs
✓ 5" Covered Castors	4pcs	✓ PP Mattress Retainer	2pcs
✓ Foldable Side Rail	1 set	☐ Basket	
✓ PP Bed Ends with Safe Lock	1 set		
✓ IV Pole 4 Hooks	1 set		

K3k

Manual Bed SK-CD1



Technical Parameters

External Length	2160mm
External Width	980mm
Bed Height Min.	400mm
Bed Height Max.	710mm
Mattress Platform	1930×830mm
Safe Working Load	200kg


Technical Configuration

✓ Manual Crank	3pcs	✓ Drainage Hook	2pcs
✓ 5" Covered Castors	4pcs	☐ Basket	
✓ Foldable Side Rail	1 set		
✓ PP Bed Ends with Safe Lock	1 set		
✓ IV Pole 4 Hooks	1 set		
✓ IV Pole Prevision	4pcs		



K2k

Manual Bed SK-C1



Technical Parameters

External Length	2160mm
External Width	980mm
Bed Height	500mm
Mattress Platform	1930×830mm
Safe Working Load	200kg

Technical Configuration

✓ Manual Crank	2pcs	✓ IV Pole Prevision	4pcs
✓ 5" Covered Castors	4pcs	✓ Drainage Hook	2pcs
✓ Foldable Side Rail	1 set	☐ Basket	
✓ PP Bed Ends with Safe Lock	1 set		
✓ IV Pole 4 Hooks	1 set		






SK-HD03

Rehabilitation Bed

320_{kg}

Maximum patient weight

The extra stable and sturdy frame allows a safe load of up to 320 kg.



SK-HD01

Electric Homecare Bed



SK-HD02

Electric Homecare Bed



SK-HD04

Electric Homecare Bed



SK-HD05

Electric Homecare Bed



Model	SK-HD01	SK-HD02	SK-HD03	SK-HD04	SK-HD05
External size	2140×1055mm	2120×1025mm	2195×1345mm	2060×1050mm	2060×1070mm
Mattress platform	2000×900mm	2000×900mm	2010×1200mm	2000×900mm	2000×900mm
Back-rest adjustment	0-72°	0-72°	0-70°	0-68°	0-68°
Leg-rest adjustment	0-28°	0-28°	0-24°	0-28°	0-28°
Reverse trendelenburg	0-12°	0-12°	0-12°	0-15°	0-17°
Trendelenburg	0-12°	0-12°	0-12°	0-18°	0-15°
Lowest & highest position	190-610mm	390-810mm	230-660mm	260-760mm	300-800mm
Safe working load	215kg	215kg	320kg	215kg	215kg
Mattress size	L2000×900mm/800mm	L2000×900mm/800mm	L2000×1200mm	L2000×900mm/800mm	L2000×900mm/800mm
Guardrail	Wooden full-length	Wooden full-length	Wooden full-length	Matel collapsible	Wooden split
Construction	Foldable	Foldable	Foldable	Non-foldable	Non-foldable
Feature	Super low position	Economy model	Super load bearing Super wide width	Enormous height range	Enormous height range

Q5k
Electric Children Bed
SKD-A



360° View


Technical Parameters

External Length	1930mm
External Width	880mm
External Height	600mm
Mattress Platform	1710×780mm

Technical Configuration

✓ Linear Motor	2pcs	✓ 5" Double Side Castors	4pcs
✓ Four-fold Aluminum Alloy Side Rails	1set	✓ IV Pole Prevision	4pcs
✓ Bed End With Medical Record Card Holder	1set	✓ IV Pole	1set
✓ Detachable PP Bed Surface With Mattress Holder	1set	✓ Drainage Hook	2pcs
✓ Hand Controller	1pc		


Q2w
Children Bed
SK-C1



Technical Parameters

External Length	2185mm
External Width	880mm
External Height	540mm
Mattress Platform	1940×780mm

CQ2k
Children Bed
SK-C2



Technical Parameters

External Length	1975mm
External Width	880mm
External Height	540mm
Mattress Platform	1730×780mm



Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole Prevision	4pcs
✓ Six-fold Aluminum Alloy Side	1set	✓ IV Pole	1set
✓ 5" Double Side Castors, Central Locking	4pcs	✓ Drainage Hook	2pcs
✓ Detachable PP Bed Surface With Mattress Holder	1set	✓ Shoes Holder	1set
✓ Bed End With Medical Record Card Holder	1set		



Technical Configuration

✓ Manual Cranks	2pcs	✓ 5" Covered Castors	1set
✓ Bed End With Medical Record Card Holder	4pcs	✓ IV Pole Prevision	4pcs
✓ Twelve-fold Foldable Stainless Steel Aluminum Alloy Side Rails	1set	✓ IV Pole	1set
✓ Detachable PP Bed Board With Built-in Mattress Holder	1set	✓ Drainage Hook	2pcs

CQ8k
Electric Children Bed
SKD-C



360° View


Technical Parameters

External Length	1975mm
External Width	880mm
Bed Height Min.	480mm
Bed Height Max.	770mm
Mattress Platform	1730×780mm

Technical Configuration

✓ LINAK Motor	4pcs	✓ PP Bed Ends with Safe Lock	1 set
✓ Hand Controller	1 pc	✓ IV Pole	1 set
✓ PP Platform with Integrated Mattress Retainer	1 set	✓ IV Pole Prevision	4pcs
✓ 5" Double Side Castors, Central Locking	4pcs	✓ Drainage Hook	2pcs
✓ Aluminium Alloy Side Rail 12 Columns with Extra Protection	1 set	✓ Electric CPR	

SK-CHD01
Super Low Children Bed



Technical Parameters


External Length	2040mm
External Width	930mm
Bed Height Min.	268mm
Bed Height Max.	688mm
Mattress Platform	1900×780mm

Technical Configuration

✓ Linear Motor	4pcs
✓ Wooden Side Rails	1 set
✓ Head & Foot Board	1 set
✓ Castors	4pcs
✓ Remote Controller	1 pc



CD3q
Children Bed
SK-CD2



Technical Parameters

External Length	2030mm
External Width	1045mm
Bed Height Min.	500mm
Bed Height Max.	770mm
Mattress Platform	1705×900mm

Technical Configuration

✓ Manual Cranks	3pcs	✓ IV Pole	1 pc
✓ 5" Hardshell Castors	4pcs	✓ Drainage Hook	2pcs
✓ Aluminium Alloy Full Length Side Rail	1 set		
✓ ABS Grace Bed Ends with Safe Lock	1 set		
✓ IV Pole Prevision	4pcs		



CX2x
Children Bed
SK-C2



Technical Parameters

External Length	1910mm
External Width	930mm
Bed Height	580mm
External Height	1260mm
Mattress Platform	1775×710mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole	1pc
✓ 5" Double Side Castors	4pcs	✓ Buffer Air Spring	4pcs
✓ Full Length Side Rail	1set		
✓ Plastic Bed Ends	1set		
✓ IV Pole Prevision	4pcs		

CR0q
Children Bed
SK-A2




Technical Parameters

External Length	1600mm
External Width	750mm
Bed Height	500mm
External Height	1040mm
Mattress Platform	1580×660mm

Technical Configuration

✓ 5" Covered Castors	4pcs	✓ Drainage Hook	2pcs
✓ Full Length Side Rail	1set		
✓ Integrated Bed Ends	1set		
✓ IV Pole	1pc		
✓ IV Pole Prevision	4pcs		

CR2q
Children Bed
SK-C2



Technical Parameters

External Length	1800mm
External Width	860mm
Bed Height	600mm
External Height	1110mm
Mattress Platform	1600×680mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole Prevision	4pcs
✓ Full Length Side Rail	1set	✓ Drainage Hook	2pcs
✓ Intergrated Bed Ends	1set	✓ Shoes Holder	1pc
✓ Anti-slip Bed Feet Cover	4pcs		
✓ IV Pole	1pc		

X01-1

Baby Crib
SK-AF3



Technical Parameters

External Length	880mm
External Width	520mm
Bed Height Min.	790mm
Bed Height Max.	960mm

Technical Configuration

✓ Gas Spring	2pcs
✓ ABS Bassinet	1pc
✓ Φ 100mm Silent Castors (2 with brakes)	4pcs

X01-5

Baby Crib
SK-AE3



Technical Parameters

External Length	810mm
External Width	520mm
Bed Height	990mm

Technical Configuration

✓ Bassinet	1pc	✓ IV Pole	1set	✓ 2cm Thickness Mattress	1pc
✓ 3" Covered Castors	4pcs	✓ Utility Basket	1pc		
✓ IV Pole Prevision	2pcs	✓ Adjustable Air Spring	1pc		

X01

Baby Crib
SK-AE3



Technical Parameters

External Length	810mm
External Width	490mm
Bed Height	960mm

Technical Configuration

✓ Bassinet	1pc	✓ IV Pole	1set
✓ 3" Castor With Brake	4pcs	✓ 2cm Thickness Mattress	1pc
✓ IV Pole Prevision	1pc	✓ Stainless Steel Frame	

X20

Pediatric Exam Table
SKZ-A

Technical Parameters

External Length	1800mm
External Width	650mm
Bed Height	830mm

Technical Configuration

✓ Material: laminate Wood Structure With PU Mattress
✓ Foot Part Can Be Slided Out For Tall Childred
✓ With Lockers At Table Side
✓ With Cartoon Stickers



SK-PB025

3D Swing Spinal Therapy Bed



Technical Parameters

External Length	1700mm
External Width	560mm
Bed Height Min.	570mm
Bed Height Max.	890mm
Safe Working Load	200kg

Technical Configuration

✓ Adjustable width headrests	2pcs	✓ Tension controlled drop devices	4pcs	✓ High-resistance springs	2pcs
✓ Fixed armrests	2pcs	✓ Stainless steel handles	4pcs	✓ Stainless steel operating lever	1pc
✓ Tissue holder	1set	✓ Bakelite adjustment knob	1set	✓ Leather bed surface	1set
✓ Electric actuator	3pcs	✓ Large hydraulic gas rod	1pc	✓ Square tube bed frame	1set
✓ Foot controller	1set	✓ Small hydraulic gas rods	2pcs	✓ Foot restraint strap	1set

SK-PB023

Electric Spinal Therapy Bed



Technical Parameters

External Length	1800mm
External Width	550mm
Bed Height Min.	550mm
Bed Height Max.	750mm
Safe Working Load	200kg

Technical Configuration

✓ Adjustable width headrests	2pcs	✓ Stainless steel handles	4pcs
✓ Fixed armrests	2pcs	✓ Bakelite adjustment knob	1set
✓ Tissue holder	1pc	✓ Leather bed surface	1set
✓ Electric actuator	3pcs	✓ Square tube bed frame	1set
✓ Foot controller	1set	✓ 3-inch castors	2pcs
✓ Tension controlled drop devices	4pcs	✓ Anti-slip foot pads	2pcs

SK-PB021

Electric Spinal Therapy Bed



Technical Parameters

External Length	1780mm
External Width	550mm
Bed Height Min.	560mm
Bed Height Max.	780mm
Safe Working Load	200kg

Technical Configuration

✓ Adjustable width headrests	2pcs	✓ Square pillow	1pc
✓ Fixed armrests	2pcs	✓ Electric actuator	2pcs
✓ Hydraulic gas rod	2pcs	✓ Foot controller	1set
✓ Tension controlled drop devices	5pcs	✓ PVC leather bed surface	1set
✓ Foot pedals	8pcs	✓ Square tube bed frame	1set
✓ Bakelite adjustment knob	1set	✓ Anti-slip foot pads	4pcs

SK-PB018

Spinal Therapy Bed



Technical Parameters


External Length	1800mm
External Width	550mm
Bed Height	530mm
Safe Working Load	200kg

Technical Configuration

✓ Adjustable width headrests	2pcs	✓ Bakelite adjustment knob	1set
✓ Fixed armrests	2pcs	✓ Durable PVC leather bed surface	1set
✓ Hydraulic gas rod	1pc	✓ Low-carbon square tube bed frame	1set
✓ Tension controlled drop devices	3pcs	✓ 3-inch wear-resistant nylon casters	2pcs
✓ Foot pedals	6pcs	✓ Anti-slip foot pads	4pcs

X36-1

Electric Examination Table



Technical Parameters


External Length	2010mm
External Width	690mm
Table Height Min.	480mm
Table Height Max.	910mm

Technical Configuration

✓ Back Gas Spring Adjustment	1pc	✓ Hand Controller	1pc
✓ Leg Gas Spring Adjustment	1pc		
✓ Domestic Single Motor	2pcs		
✓ 3-inch Universal Wheels	4pcs		

X36

Electric Examination Table



Technical Parameters

External Length	1860mm
External Width	680mm
Table Height Min.	490mm
Table Height Max.	910mm

Technical Configuration

✓ Back Gas Spring Adjustment	1pc	✓ Hand Controller	1pc
✓ Leg Gas Spring Adjustment	1pc	✓ Pillow	1pc
✓ Domestic Single Motor	1pc		
✓ 3-inch Universal Wheels	4pcs		

X35

Electric Examination Table



Technical Parameters


External Length	1845mm
External Width	680mm
Table Height Min.	550mm
Table Height Max.	860mm

Technical Configuration

✓ Backrest Lifting Rod Adjustment	1pc	✓ Paper Roll Holder	1pc
✓ Domestic Single Motor	1pc	✓ Hand Controller	1pc
✓ 2-inch Hidden Castors	2pcs	✓ Pillow	1pc
✓ 2-inch Universal Wheel	2pcs		

X14

Hydraulic Examination Table SKZ-BD



Technical Parameters

External Length	1940mm
External Width	700mm
Table Height Min.	490mm
Table Height Max.	890mm

Technical Configuration

✓ Hydraulic Pump	1pc	✓ Toilet Paper Holder	1pc
✓ Adjustable Air Spring	2pcs		
✓ 5" Covered Castors	4pcs		
✓ Hydraulic pedal	2sets		

X07a
Examination Table
SKZ-A



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	680mm

Technical Configuration

✓ Anti-slip Bed Feet Cover	4pcs

X09
Examination Table
SKZ-B




Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	680mm

Technical Configuration

✓ Anti-slip Bed Feet Cover	4pcs
✓ Mechanical Lever to Adjust Back Rest	1set

A048
Gynecological Examination Table
SKZ-B



Technical Parameters

External Length	1810mm
External Width	690mm
Table Height	880mm

Technical Configuration

✓ Drawer	3pcs	✓ Lithotomy Rod	1set
✓ Step Stool	1pc	✓ Stainless Steel Basin	1pc
✓ Retractable Foot Board	1pc	✓ Paper Roll Holder	1pc

A98-8
Electric Obstetric Bed




Technical Parameters

External Length	2300mm
External Width	1050mm
Table Height Min.	480mm
Table Height Max.	880mm

Technical Configuration

✓ Guardrail	2pcs	✓ Delivery Holder	1pc	✓ Mattress	1set
✓ Head&foot Board	1pc	✓ Foot Support	1pc	✓ Basin	1pc
✓ Motor	1set	✓ 6" Castor, Central Locking System	4pcs	✓ IV Pole	1pc
✓ Obstetric Handle	2pcs	✓ Wall Bumper Wheel	2pcs	✓ Backup Battery	1pc
✓ Auxiliary Table	1pc	✓ Night Light	1pc		

A98-3
Electric Obstetric Bed



Technical Parameters

External Length	2120mm
External Width	1120mm
Table Height Min.	595mm
Table Height Max.	905mm

Technical Configuration

✓ Pull Handle	2pcs	✓ Stainless Steel Basin	1pc
✓ Leg support	2pcs	✓ Instrution Manual	1pc
✓ Foot support	2pcs	✓ 6" Castor	1pc
✓ Hand Remote Controller	1pc	✓ Mattress	1set

A98-4

Electric Obstetric Bed

360mm
Height
Adjustment

70°
Back-rest

10°
Reverse
Trendelenburg

10°
Trendelenburg

65°/40°
Leg Plate
up/down



Technical Parameters

External Length	2050mm
External Width	960mm
Table Height Min.	640mm
Table Height Max.	1000mm

Technical Configuration

✓ Linear Motor	3pcs	✓ Central Locking Pedal	2pcs	✓ Hand Controller	1pc
✓ 5" Double Side Castors, Central Locking	4pcs	✓ Hand Grip	1 pair	✓ Power Line	1pc
✓ ABS Side Rail With Angle Indicator	1set	✓ Leg Holder	1 pair	✓ Leg Strip	1 Pair
✓ ABS Head Board With Safe Lock	1pc	✓ Basin	1pc	✓ IV Pole Prevision	2pcs

A99-13

Electric Obstetric Bed

250mm
Height
Adjustment

50°
Back Plate
up

22°
Reverse
Trendelenburg

22°
Trendelenburg

90°
Leg Plate
Swing Out

16°
Seat Plate
up



Technical Parameters

External Length	1860mm
External Width	600mm
Table Height Min.	730mm
Table Height Max.	980mm

Technical Configuration

✓ Armrest	2pcs	✓ SS Rinsing Basin	1pc
✓ Knee Curtch	2pcs	✓ Pedal Controller	1pc
✓ Pull Handle	2pcs	✓ LED Examination Light	1pc
✓ Anesthesia Frame	1pc		

A99-7

Electric Obstetric Bed

250mm
Height
Adjustment

55°/18°
Back Plate
up/down

35°/3°
Seat Plate
up/down

90°
Leg Plate
Swing Out



Technical Parameters

External Length	1980mm
External Width	720mm
Table Height Min.	630mm
Table Height Max.	880mm

Technical Configuration

✓ Bed Frame	1pc	✓ Hand Bars	2pcs	✓ Castor	1pc
✓ Foot Pedal	1pc	✓ Leg Support with Cushion	2pcs		
✓ Hand Controller	1pc	✓ Arm Support with Cushion	2pcs		
✓ Anesthesia Screen	1pc	✓ Stainless Steel Basin	1pc		

A99-5
Electric Obstetric Bed

300mm
Height
Adjustment

55°/10°
Back Plate
up/down

20°/10°
Seat Plate
up/down

90°
Leg Plate
Swing Out



Technical Parameters

External Length	1985mm
External Width	665mm
Table Height Min.	675mm
Table Height Max.	975mm

Technical Configuration

✓ Arm Supports with Cushions	2pcs	✓ Hand Controller	1pc
✓ Hand Grips (Side Rails)	2pcs	✓ Foot Controller (Pedal Switch)	1pc
✓ Leg Supports with Cushions	2pcs	✓ Fixed mattress 65mm	1pc
✓ Stainless Steel Waste Basin	1pc	✓ Built-in battery	1pc

A105
Manual Obstetric Table

250mm
Height
Adjustment

50°/3°
Back Plate
up/down

10°
Reverse
Trendelenburg

20°
Trendelenburg

90°
Leg Plate
Swing Out



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	700mm
Table Height Max.	950mm

Technical Configuration

✓ Bed Frame	1pc	✓ Arm Support with Cushion	2pcs
✓ 50mm Mattress	1pc	✓ Stainless Steel Basin	1pc
✓ Hand Bars	2pcs	✓ Auxiliary plate	1pc
✓ Leg Support with Cushion	2pcs		

A99-8
Electric Gynecological Exam Couch

70°/20°
Back Plate
up/down

40°
Seat Plate
up

490mm
Height
Adjustment



Technical Parameters

External Length	1400mm
External Width	740mm
Table Height Min.	610mm
Table Height Max.	1100mm

Technical Configuration

✓ Foot Switch	1pc	✓ Removable Handles	2pcs
✓ Power Switch	1pc	✓ Pull-out Stainless Steel Basin	1pc
✓ Seamless Foaming Mattress	1set	✓ 5" Double Sided Castor	1set
✓ Leg Support	2pcs		

A045-1
Gynecological Table
SKZ-C

70°±10°
Back Plate
up

90°±10°
Leg Plate
down



Technical Parameters

External Length	1820mm
External Width	590mm
Table Height	800mm

Technical Configuration

✓ Leg Support	2pcs	✓ Mechanical Lever	2sets
✓ Hand Bar	1set	✓ Stainless Steel Basin	1pc
✓ Anti-slip Bed Feet Cover	4pcs		



A307

Electric-Hydraulic Operating Table

The A307 operating table is made of high-strength synthetic plate which is translucent for X-ray. It can also shift the table back and forth up to 300mm for C-arm examination of all parts of the patient's body during the operation without moving the patient.

30°
Reverse
Trendelenburg

30°
Trendelenburg

20°
Lateral Tilt
left/right

45°/90°
Head Plate
up/down

80°/50°
Back Plate
up/down

15°/90°
Leg Plate
up/down

90°
Leg Plate
Horizontal
Positioning

300mm
Horizontal
Sliding
C type Arm

100°
Flex

220°
Reflex



Technical Parameters

External Length	2140mm
External Width	550mm
Table Height Min.	550mm
Table Height Max.	1050mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair		
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		



A302
Electric Operating Table
SKL-A

A302 is a highly automated electric Multi-functional operating table. Its elevation, lowering, trendelenburg, reverse trendelenburg, lateral tilt, back folding, vertical translation and leg folding are all electrically operated. The operating table is equipped with a digital display hand controller to make operation easy and intuitive. With Imported large diameter castors, smooth movement, stable brake system. The standard dual control system guarantees the safety of operation.

24°
Reverse
Trendelenburg

24°
Trendelenburg

23°
Lateral Tilt
left/right

33°/90°
Head Plate
up/down

68°/22°
Back Plate
up/down

18°/85°
Leg Plate
up/down

110mm
Kidney
Bridge Up

350mm
Horizontal
Sliding
C type Arm

100°
Flex

220°
Reflex



Technical Parameters

External Length	2000mm
External Width	555mm
Table Height Min.	690mm
Table Height Max.	990mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		



A301
Electric Operating Table
SKL-A

A301 is a multi-functional electric operating table. Its mattress is made of polyurethane molded material, can be sterilized and cleaned at high temperature, and there are no gaps so to avoid bacterial growth. The platform is made of high-purity phenolic resin board applicable for C-arm, and can be slid up to 350mm in one direction. The X-ray photography can carry out the whole body radiography examination without moving the patient, no blind spots or shadows.

24°
Reverse
Trendelenburg

24°
Trendelenburg

23°
Lateral Tilt
left/right

30°/90°
Head Plate
up/down

66°/22°
Back Plate
up/down

90°
Leg Plate
down

90°
Leg Plate
Horizontal
Positioning

110mm
Kidney
Bridge Up

350mm
Horizontal
Sliding
C type Arm



Operating Table

Technical Parameters

External Length	2000mm
External Width	555mm
Table Height Min.	730mm
Table Height Max.	980mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		

A100-4A

Electric Operating Table

SKL-A

A100-4A is an electric operating table with super load-bearing and extra-long sliding distance , which can be applicable for C-arm. The built-in large-capacity battery can continue to support all operations in the event of a power failure. And the operating table adopts a dual-control system to ensure the safety of the operation.



Technical Parameters

External Length	2115mm
External Width	550mm
Table Height Min.	710mm
Table Height Max.	1010mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports (with straps)	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests (with straps)	1 pair		
☒ Lumbar Supports	1 pair		

A100-4

Electric Operation Table

SKL-A

A100-4 is an electric operating table that can meet the needs of various surgeries. It is driven by low-voltage DC motor, with safe and stable performance. Super wide table top, super long sliding distance. It can be applicable for C-arm. It has a high level of automation, low noise, strong reliability, and imported original parts for key parts. It is a perfect and ideal operating table.



Technical Parameters

External Length	2080mm
External Width	550mm
Table Height Min.	700mm
Table Height Max.	1000mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports (with straps)	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests (with straps)	1 pair		
☒ Lumbar Supports	1 pair		

Operating Table

A501

Electric-Hydraulic Operating Table

A501 is an advanced, multi-functional electric-hydraulic operating table. Its table top is made of high-strength synthetic board and can transmit X-rays. In addition, the platform can be equipped with a 180° rotation function to realize head and leg interchange, which is applicable for C-arm perspective and placement of some special positions.



Technical Parameters

External Length	2125mm
External Width	550mm
Table Height Min.	630mm
Table Height Max.	980mm


Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair		
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		

A2000

Electric Operating Table SKL-A

A2000 is a basic electric operating table with a very high cost effective, it can fulfill the requirement of general surgeries. Its built-in battery can continue to support all operations in the event of a power failure.



Technical Parameters

External Length	2050mm
External Width	480mm
Table Height Min.	680mm
Table Height Max.	930mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests	1 pair		
☒ Lumbar Supports	1 pair		

AM052C-1

Manual Operating Table SK1F-R



Technical Parameters

External Length	2040mm
External Width	480mm
Table Height Min.	790mm
Table Height Max.	1030mm

Technical Configuration

☒ Anesthesia Screen	1 unit	☒ Leg Supports (with straps)	1 pair
☒ Shoulder Holder	1 pair	☒ Foot Plates	1 pair
☒ Arm Rests (with straps)	1 pair		
☒ Lumbar Supports	1 pair		

Multi-functional Orthopedic Frame



Gel Mattress



Carbon Fiber Multi-functional Orthopedic Frame



Upper Limb Traction Device



Multi-purpose Head Rack



Shoulder Traction Frame



Stirrup Pneumatic Power Leg Rack



Bow Spine Operating Supporter



Mechanical Power IV Pole



Operating Table Accessories Cart



Knee Surgery Rack



Double Table Top



Prone Body Supporting Frame



Upper Limb Operating Table



Bed Surface Widening Module



External Kidney Bridge



Dirt Basin



Body Safe Strap





SHADOW MANAGEMENT


Everything has been designed to enable the medical team to focus exclusively on the surgical operation thanks to the perfect shadow dilution. A large number of individual divergent light beams create a complete light patch, ensuring that the light field will remain circular and shadow-free. Alno Light20 provides the best illumination, even when surgeons are blocking the LED.

Alno LED
operation lamp



Bulb

OSRAM LED bulbs with 72 beads and 24 triangular lenses per head, providing long-lasting, bright, and uniform illumination.



Color Temperature&illumination

Advanced LED cold light with up to 160,000 Lux, 10-level adjustment, and 3500–5500K color temperature for soft, non-glare illumination.



5-inch Extra-large LCD

5-inch large LCD with simple interface, 30° adjustable angle, and seamless design for easy disinfection.

OSRAM ^{Hours} 60000

^{Kelvin} 3500-5500

^{Angle Adjustment} ±30°

The head of the lamp is windmill-shaped, Saikang's unique (appearance patent) design. The back is equipped with a soft ambient light. The top offers a soft ambient light, simple and stylish while also being useful for low-light surgery (laparoscopic). Exceeding the infection control standards of modern operating rooms.



User friendly for multiple surgeons and nurses at the same time



Head insert design: No heating over the surgeons head (max. Allow is 1C°)

Alno /series
LED Operation Lamp



Operating Handle

Detachable PPSU operating handle, heat- and pressure-resistant for safe steam sterilization.

Intelligent Shadow Management

Smart sensors adjust lighting to eliminate shadows and maintain clear visibility.

Lighting Mode

It has three lighting modes (standard, endoscopic, and low-light), as well as red light compensation, green light compensation, and shadow compensation.



STANDARD

ENDOSCOPIC

LOW-LIGHT

Alno Light20



Alno Light10



Alno Light M




	Alno Light20 SK-SD700/700	Alno Light10 SK-SD700	Alno Light M
Central illuminance (Ec)	40000-160000 lx (10 levels)	40000-160000 lx (10 levels)	40000-160000 lx (10 levels)
Spot diameter (d10)	200-300mm (5 levels)	200-300mm (5 levels)	200-300mm (5 levels)
Spot distribution (d50/d10)	55%	55%	55%
Depth of the light column (60% of Ec)	700mm	700mm	700mm
Color temperature (Tc)	3500-5500K (11 levels)	3500-5500K (11 levels)	3500-5500K (11 levels)
Colour rendering index (Ra)	98	98	98
Maximum total irradiance (Ee)	≤550 W/m²	≤550 W/m²	≤550 W/m²
Radiant energy (Ee/Ec)	≤3.5mW/(m².lx)	≤3.5mW/(m².lx)	≤3.5mW/(m².lx)



Alno X71

Medical Bridge Ceiling Pendant (separate wet and dry areas)

Modern ceiling pendant with separate wet & dry areas, offering flexible platforms, drawers, and infusion holders for efficient and safe medical rescue.



Zoned Lighting System

Smart zoned ICU lighting with adjustable brightness and color; all controlled by touch panel.



Information Display

Built-in bed info display for quick patient data access and improved efficiency.

Alnno / series

Suitable for ICU, NICU, PICU, EICU room

Wet and Dry Areas Separated

Separated wet and dry areas ensure clear layout, minimize interference, and enhance safety.



Translation System

Translation system with 500mm range, smooth control, and precise positioning for multi-directional care.



Suspension System

Both the tower and arms can rotate 340°, with built-in limit devices to prevent drift.

Rotate
340°



Telescopic Extension Arms

Telescopic extension arms with high load capacity, flexibly adapt to infusion needs and optimize space.



Braking System

One-button remote braking with indicator lights ensures quick, precise control.



Drawer

Durable 30kg drawer with smooth silent rails and optional sensor light for convenience.



Cable Management



Cable Hub



Stainless steel basket



Suction tube box (Optional)



Drawer divider box (Optional)



Computer keyboard drawer (Optional)



Flashlight + flashlight placement box (Optional)



Stethoscope + stethoscope hanger (Optional)

Alnno L21

Electric Double-Arm Pendant

- Lightweight, enclosed design – safe and flexible.
- Electric double-arm – wide treatment coverage.
- Supports multiple devices – space-efficient.
- Boosts OR efficiency – saves space.




Alnno X22

Double-Arm Pendant

- Lightweight, compact design – easy and safe to use.
- Double-arm movement – extends treatment reach.
- Accommodates multiple devices – maximizes workspace.
- Enhances OR workflow – saves space and time.




Mattress (OEM Size)



SKP028

Pressure Relieving Mattresses with Gel




SKP040

Cube Design Pressure-relieving Mattress



SKP041

Pressure-relieving Mattress



Mattress (OEM Size)

SKP011




High Density Foam Memory Foam



PU Composite



SKP001



Air Through Water Proof Sectional Foam



PU Composite



TPU Fabric Composite



Canvas Composite



SKP001



SKP002



SKP003

SKS009-2

Phenolic Resin Bedside Table

Technical Parameters

External Size	590x480x840mm
Dining Table Height Adjustment	750-1000mm

Technical Configuration

✓ Top Guardrail	Aluminum alloy
✓ Drawer	1pc
✓ Cabinet	1pc
✓ Casters	4pcs
✓ Side Pull-up Writing Board	1pc



SKS02-W

ABS Bedside Table

475×465×795mm(L×W×H)



SKS03-W

ABS Bedside Table

480×450×950mm(L×W×H)



SKS007-4

Steel Bedside Table

440×460×810mm(L×W×H)



SKS036

Bedside Table

440×440×790mm(L×W×H)



SKS008

ABS+Steel Bedside Table

500×510×705mm(L×W×H)



SKS013-1

ABS+Steel Bedside Table

575×500×750mm(L×W×H)





SKS025
Phenolic Resin Bedside Table
475×475×800mm(L×W×H)

SKH046-11
Overbed Table

- Technical Parameters
800×395×870-1170mm
- Technical Configuration
- ✓ Air Spring for Height Adjustment
 - ✓ Phenolic Resin Table Top



SKH046-2
Overbed Table



- Technical Parameters
890×445×760-1050mm
- Technical Configuration
- ✓ Air Spring for Height Adjustment
 - ✓ PE Table Top with Cup Socket



SKH046-14
Overbed Table



- Technical Parameters
890×400×890-1190mm
- Technical Configuration
- ✓ Air Spring for Height Adjustment
 - ✓ ABS Table Top with Cup Socket and Push Handle



SKH042
Overbed Table



- Technical Parameters
800×400×810-1110mm
- Technical Configuration
- ✓ Air Spring for Height Adjustment
 - ✓ MDF Table Top



SKH046-SC
Telescopic Dining Board



- Technical Parameters
890-1130×350mm
- Technical Configuration
- ✓ Hook



SKE-136
Electric Dalysis Chair
SK1D-0D2



Technical Parameters


External Length	2020mm
External Width	600mm
Chair Height Min.	620mm
Chair Height Max.	870mm

Technical Configuration

✓ Hand Controller	1pc	✓ Arm Support	2pcs
✓ USB Port	1pc	✓ Plug-in Headrest (with a Small Pillow)	1pc
✓ 3" Brake Wheel	4pcs	✓ Seamless Soft - Padded Seat Cushion	1set
✓ Small Dining Table	1pc		



SKE-132
Electric Dalysis Chair



Technical Parameters

External Length	1730mm
External Width	610mm
Min seat height from ground	520mm
Max seat height from ground	550mm

Technical Configuration

✓ Dining Tray	1set
✓ Castors	4pcs
✓ Armrests	2pcs
✓ Hand Controller	1pc



SKE001

Attendant Chair

Technical Parameters

Seat Size (L×W×H)	680×720×940mm
Bed Size (L×W×H)	680×1900×590mm



SKE001-5

Therapy Chair

Technical Parameters

Seat Size (L×W×H)	650×780×1060mm
Bed Size (L×W×H)	650×2060×450mm



SKE001-16

Attendant Sofa

Technical Parameters

Seat Size (L×W×H)	800×1000×830mm
Bed Size (L×W×H)	800×1960×500mm



SKE001-19

Attendant Sofa

Technical Parameters

Seat Size (L×W×H)	1600×700×655mm
Bed Size (L×W×H)	2200×700×655mm



SKE091

Blood Donation Chair

Technical Parameters

Seat Size (L×W)	1750×600/940mm
Seat Height	480-540mm



SKE941

Therapy Chair

Technical Parameters

Overall Size	770×940×1250mm	Seat Height	560mm
Seat Cushion Size	530×480mm	Back-rest	10°-70°



SKE001-10

Attendant Sofa

Technical Parameters

Seat Size (L×W×H)	850×1000×850mm
Bed Size (L×W×H)	850×1950×630mm



SKE014-9

Attendant Sofa

Technical Parameters

Seat Size (L×W×H)	1250×1000×850mm
Bed Size (L×W×H)	1250×1950×630mm





SKE909

Reclining Sofa

Technical Parameters

Seat Size (L×W×H)	640×900×1290mm
Bed Size (L×W)	640×1950mm



SKE087

Reclining Chair

Technical Parameters

Back size (L×W)	710×540mm
Seat Size (L×W)	450×540mm
Leg Size (L×W)	540×540mm



SKE090

Blood Donation Chair

Technical Parameters

Seat Size (L×W×H)	905×640×920mm
-------------------	---------------

Technical Configuration


- ✓ Detachable frame made of powder coated steel
- ✓ PU cover with foam
- ✓ The armrest can rotate 90 degrees and The height can be adjusted
- ✓ The nearby drawer cabinet can be used for storage
- ✓ Rubber feet pad



SKE011-1

Treat-waiting Chair


Material: PU
1630×650×865mm(L×W×H)



SKE011-6

Treat-waiting Chair


Material: PU
1655×640×820mm(L×W×H)



SKE011-8

Treat-waiting Chair

Material: PU
1780×600×820mm(L×W×H)





360° View

SK-CL001 SK-CL501
Curtain Track



SKH049-7
Phenolic Resin Screen
Full Size: 1780×1650mm
Folded Size: 594×412×175mm



SKH049
Stainless Steel Screen
Size: 2000×1870mm



DM-Z104
Office Table
Size: 1400×600-1040×750mm



DM-Y106
Office Chair
Size: 650×640×965mm



SKE801
Nurse Stool
Size: Φ360×450-570mm



SKE020-1
Step Stool
Size: 490×410×170/340mm



SKH034
Trash Can



12L

SKH034-5
Kick Bucket
Size: Φ400×H360×14L



14L

SKH041-3 I.V. Stand

Function

Height Adjustment 1350-2250mm

Material

IV Pole Outer Tube: Φ25×1000mm (Stainless Steel)

IV Pole Inner Tube: Φ16×1200mm (Stainless Steel)



SKH041(1) I.V. Stand

Function

Height Adjustment: 1350-2250mm

Material

IV Pole Outer Tube

Φ25×1000mm (Stainless Steel)

IV Pole Inner Tube

Φ16×1200mm (Stainless Steel)



SKH041(18) I.V. Stand

Function

Height Adjustment: 1200-2000mm

Material

IV Pole Outer Tube

Φ32mm (Stainless Steel)

IV Pole Inner Tube

Φ19mm (Stainless Steel)



SKH041(16) I.V. Stand

Function

Height Adjustment: 1500-2000mm

Material

IV Pole Outer Tube

Φ32mm (Stainless Steel)

IV Pole Inner Tube

Φ22-Φ25mm (Stainless Steel)



SKB1C10
Electric Stair Stretcher



Technical Parameters

Full Size	1225×510×1620mm
Folded Size	1080×510×270mm

SKB039(H)
Stretcher Trolley
SKC-3B4



Technical Parameters

Full Size	1950×570×910mm	Mattress thickness	PVC, 100 mm
Folded Size	2030×570×300mm	Load Bearing	180kg
Back-rest	75°		

SKB039(G)
Stretcher Trolley
SKC-3B10



Technical Parameters

Full Size	2000×560×1000mm	Knee-rest	20°
Folded Size	2000×560×380mm	Mattress thickness	PVC, 100 mm
Back-rest	75°	Load Bearing	200kg


SKB039(F)
Stretcher Trolley
SKC-3F5



Technical Parameters

Full Size	1920×570×825mm	Back-rest	80°
Folded Size	2000×570×260mm	Load Bearing	280kg

SKB037(B)
Patient Trolley
SKM-C



Technical Parameters

Length	1960mm
Width	520mm
Height	770mm

Technical Configuration

✓ 5 inch wheel:	4pcs
✓ Mechanical lever (10 sections)	1set
✓ Stainless steelw pole(Φ13-16)	1pc

SKB041-3S
Hydraulic Patient Trolley
SKM-E



Technical Parameters


Length	2150mm
Width	790mm
Min. Height	635mm
Max. Height	895mm

Technical Configuration

✓ Hydraulic Pump	1pc	✓ Wall bumper wheel	4pcs	✓ Debris basin	1pc
✓ Air Spring	1pc	✓ Foot pedal	2pcs	✓ IV pole	1pc
✓ Stainless steel guardrail	1set	✓ 5" Castors with brake	4pcs	✓ Mattress with thickness of 60mm	1pc
✓ Push handle	1set	✓ Oxygen cylinder rack	1pc		



SKB041-7
Patient Transportation Trolley
SKV-BE



Technical Parameters

Length	2200mm
Width	770mm
Height (X-ray Backrest)	550-860mm
Height (X-ray Platform)	600-910mm

Technical Configuration

☒ 8-inch Central Control Castors	4pcs	☒ Leg Lifting Crank Handle	1pc	☐ Paper Holder	1pc
☒ Hydraulic Pumps	2pcs	☒ Anti-collision Wheels	4pcs	☐ Foot-Board & Monitor Holder	1pc
☒ 5-speed Foldable Guardrails	1pair	☒ Foot Pedals	2sets	☐ Monitor Holder	1pc
☒ Stainless Steel Infusion Pole	1pc	☒ Central Control Brakes	4pcs	☐ Debris Basket	1pc
☒ Infusion Stand Holes	4pcs	☒ Fifth Castor	1pc	☐ Weighting System	1pc
☒ Back Lifting Gas Springs	2pcs	☒ Bottom Plastic Cover	1pc	☐ Pluggable Mobile IV pole	1set
☒ Back Angle Displays	2pcs	☒ 8cm PU Mattress	1pc		
☒ Patient Restraints/Straps	2pcs	☐ Oxygen Tanks Holder	1pc		

SKB041-1

Patient Transportation Trolley

SKM-D



Technical Parameters

Length	1900mm
Width	640mm
Min. Height	550mm
Max. Height	850mm

Technical Configuration

✓ Crank System	1 set	✓ Directional Wheel	1 pc	✓ Central Locking Break Pedal	2 pcs
✓ Back Rest Gas Spring	1 pc	✓ IV Pole Previsions	2 pcs	✓ Directional Pedal	2 pcs
✓ ABS Side Rail	1 set	✓ IV Pole	1 set	✓ 6cm PU Mattress	1 pc
✓ 6-inch Castors with Central Locking	4 pcs	✓ Oxygen Cylinder Holder	1 pc		

SKB041-2

Patient Transportation Trolley

SKV-BC



Technical Parameters

Length	1910mm
Width	640mm
Min. Height	650mm
Max. Height	960mm

Technical Configuration

✓ 8-inch Central Control Castors	4 pcs	✓ PP Guardrails	2 pcs	✓ Fifth Castor	1 pc
✓ Hydraulic Pumps	2 pcs	✓ Back Gas Spring	1 pc	✓ Bottom Plastic Cover	1 pc
✓ Stainless Steel Infusion Pole	1 pc	✓ Foot Pedals	2 sets	✓ 6cm PU Mattress	1 pc
✓ Infusion Stand Holes	2 pcs	✓ Central Control Brakes	2 pcs		

SKB041-6

Patient Transportation Trolley

SKM-F



Technical Parameters

Length	1930mm
Width	663mm
Min. Height	510mm
Max. Height	850mm

Technical Configuration

✓ Back Rest Air Spring	1 pc	✓ IV Pole Previsions	4 pcs	✓ Central Locking Break Pedal	4 pcs
✓ Lifting Screw Rod	1 pc	✓ IV Pole	1 pc	✓ Bottom Plastic Cover	1 pc
✓ PP Material Bed Boards	2 pcs	✓ Fifth Wheel Device	1 set	✓ Oxygen Cylinder Holder	1 pc
✓ PP Material Guardrails	2 pcs	✓ 6-inch Central Control Wheels	4 pcs	✓ 6cm PU Mattress	1 pc

SKB041-10

Patient Transportation Trolley

SKV-BE



Technical Parameters

Length	2150mm
Width	830mm
Min. Height	560mm
Max. Height	890mm

Technical Configuration

✓ 8-inch Central Control Castors	4 pcs	✓ Back Lifting Gas Springs	2 pcs	✓ Fifth Castor	1 pc
✓ Hydraulic Pumps	2 pcs	✓ Leg Lifting Crank Handle	1 pc	✓ Bottom Plastic Cover	1 pc
✓ PP Guardrails	2 pcs	✓ Anti-collision Wheels	4 pcs	✓ 7cm PU Mattress	1 pc
✓ Stainless Steel Infusion Pole	1 pc	✓ Foot Pedals	2 sets	□ Instrument Platform	1 pc
✓ Infusion Stand Holes	4 pcs	✓ Central Control Brakes	4 pcs		

SKR-R10

ECG Machine Cart



SKR-R03

Laptop Cart



SKR-R04

Mobile Charting Cart



SKR-R07

Mobile Charting Cart



Model	SKR-R10	SKR-R03	SKR-R04	SKR-R07
Length	500mm (±10mm)	530mm (±10mm)	535mm (±10mm)	530mm (±10mm)
Width	500mm (±10mm)	520mm (±10mm)	520mm (±10mm)	620mm (±10mm)
Height	875mm (±10mm)	725-1045mm (±10mm)	725-1045mm (±10mm)	855-1100mm (±10mm)
Tabletop Size	335×365mm (±10mm)	450×445mm (±10mm)	450×505mm (±10mm)	480×460mm (±10mm)
Auxiliary Workbench	/	✓	/	/
Keyboard Case	/	✓	/	✓
Mouse Hanger	/	/	/	✓
Storage Box	/	/	/	✓
Column Material	Aluminum Alloy	Aluminum Alloy	Aluminum Alloy	Aluminum Alloy
Power System	/	Pneumatic Lifting	Pneumatic Lifting	Pneumatic Lifting
Base	Carbon Steel Base + ABS Cover	Cast Aluminum with Spraying Coating	Cast Aluminum with Spraying Coating	Cast Aluminum with Spraying Coating
Castor	5×2.5-inch double-sided castors, 2 of which are with brakes	4×4-inch double-sided castors, 2 of which are with brakes	4×4-inch double-sided castors, 2 of which are with brakes	4×4-inch double-sided castors, 2 of which are with brakes
Battery	/	/	/	Optional



SK-ET750

Emergency Trolley
Size: 750×530×1050mm

Technical Configuration

✓ Defibrillator shelf	1pc	✓ Castor	4pcs	✓ CPR Board	1pc
✓ Sliding Side Shelf	1pc	✓ Trash Can	2pcs	✓ Oxygen Bottle Holder	1pc
✓ Dust Baskets	1pc	✓ Central Lock	1pc		
✓ IV pole	2pc	✓ External Connectors	1 set		



SKR054-ET
Emergency Trolley
Size: 760×475×945mm



SKR-ET625
Emergency Trolley
Size: 625×475×945mm

SKR054-AT
Anesthesia Trolley
Size: 760×475×945mm



SKR054-MT
Medicine Trolley
Size: 750×475×920mm



SKR054-ET77A
Emergency Trolley
Size: 750×475×950mm



SKR054-AT77A
Emergency Trolley
Size: 750×475×950mm



SKR-ET625-1
Anesthesia Trolley
Size: 620×475×945mm



SKR-AT625
Anesthesia Trolley
Size: 620×475×945mm



SKR-AT625-1
Anesthesia Trolley
Size: 620×475×945mm



SKR-MT625
Medicine Trolley
Size: 620×475×945mm



SKR058-ET

Emergency Trolley
Size: 790×600×1000/1100mm



SKR058-MT

Anesthesia Trolley
Size: 790×600×940/1020mm



SKR058-NT

Nursing Trolley
Size: 620×600×1000/1100mm



SKH040

Laundry Trolley
Size: 800×600×840mm



SKH040-1G

Laundry Trolley
Size: 800×600×865mm



SKH027

Laundry Trolley
Size: 900×450×1050mm



SKH006C

Instrument Trolley
[Disassembly structure]
Size: 670×450×900mm



SKH001

Treatment Trolley
Size: 670×450×900mm




SKH004

Instrument Trolley
Size: 670×450×900mm




SKH005

Fan-shaped Instrument Trolley
Size: L: 1550×460×900mm




SKH036


Stainless Steel Washing Sink
Size: 1200(W)×600(D)×1100(H)mm



Size: 600(W)×600(D)
× 1100(H)mm



Size: 1800(W)×600(D)
× 1100(H)mm



Size: 2400(W)×600(D)
× 1100(H)mm



Technical Configuration

✓ Defibrillator shelf	1pc	✓ IV pole	1pc	✓ Sealed Handle	1pc
✓ Sliding Side Shelf	1pc	✓ Castor	4pcs		
✓ Writing Board	1pc	✓ Oxygen Bottle Holder	1pc		
✓ CPR Board	1pc	✓ Trash Can	1pc		

NOTE