

Product Datasheet

anti- FGFR1 antibody (MDB1100638)

Description	<i>anti- FGFR1 antibody...</i>
Species/Host	<i>Mouse</i>
Reactivity	<i>Human, Mouse, Rat</i>
Tested Applications	<i>ELISA, WB, IHC, IHC-P</i>
Immunogen	<i>fibroblast growth factor receptor 1</i>
Form/Appearance	<i>liquid</i>
Storage	<i>PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20 °C for 12 months(Avoid repeated freeze / thaw cycles.)</i>
Note	<i>For research use only.</i>
Isotype	<i>IgG</i>
Clonality	<i>Monoclonal</i>
Purity	<i>≥95% as determined by SDS-PAGE</i>
MW	<i>45-55 kDa</i>
Uniprot ID	P11362
NCBI	468
Dilution Range	<i>WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IHC(P): 1:50 – 1:100</i>
Expiration Date	<i>12 months from date of receipt.</i>

If you have any question on the order & technical assistance please contact us via:

Jl. Bandara Abdul Rahman Saleh GG 3-4 Malang, INA.

Website: www.medikbioindo.com

Tel: (62) 341 3021374, Fax: (62) 341 3021374 · Email: info@gbi-medika.com. All rights reserved